



ATP TESTING KIT

SystemSURE Plus ATP Hygiene Monitoring System

ATP testing is the best method of hygiene monitoring, providing results in 15 seconds of inserting the sample into the device. The use of small pen sized sample swab makes taking the sample as easier and time efficient as ever.

ATP stands for Adenosine Triphosphate and is the energy molecule found in all living cells. Testing for this indicates the amount of living cells on a surface and therefore can determine if a surface is clean or not. A special liquid is introduced into the sample, and then causes the ATP to emit light. This light is measured and is proportional to the amount of living matter taken within the sample.

TYPICAL USES:

- Hospital Cleaning- Cleanliness Indication**
 Hospitals require rigorous cleanliness standards, and in some instances need to be completely sterile. ATP testing can be a great indicator to the cleanliness of a surface post hospital cleans.
- Sewage Cleaning- Cleanliness Indication**
 Sewage cleaning requires a proper clean to ensure the safety of everyone, using ATP testing can be a great indicator for the effectiveness of the cleaning process
- Food Hygiene- Cleanliness Indication**
 A typical use for measuring cleanliness in the Hospitality and food industry, this can be a great way to indicate whether your kitchen is up to cleanliness standard
- Cleaning Process Effectiveness Indicator**
 An instantaneous way to verify if your cleaning processes are up to standard

Specifications

Power Supply	2 x AA Batteries
Weight	258g
Consumables	Single Use Sample Swabs
Dimensions	7 x 18 x 3 cm

Note

The ATP Testing Kit will consist of the following:

- 1 x System SURE Plus ATP Tester
- 2 x Sample Swabs (More can be purchased from our products page)